

# Seraseq<sup>®</sup> ctDNA v4 Reference Materials

Accelerate liquid biopsy assay validation and clinical implementation with a comprehensive reference material

## INTRODUCTION

As liquid biopsy assays gain acceptance for a growing variety of clinical applications, an acute need for highly multiplexed reference materials bearing a broad range of clinically relevant variant types has arisen. But searching for useful remnant patient specimens is prohibitively time-consuming, and singleplex cell lines are inefficient and add unnecessary expense.

Seraseq ctDNA v4 reference materials are purpose-built to accelerate the thorough validation and clinical implementation of NGS-based ctDNA assays. Seraseq ctDNA v4 contains an unprecedented 93 variants including SNPs, indels, translocations and CNVs. Seraseq ctDNA v4 is manufactured with a new fragmentation process which results lower background than previous Seraseq ctDNA offerings to support greater confidence at low variant allele frequencies. These orthogonally validated truth-set materials can trim weeks or even months off your assay validation timeline with the most complete coverage of actionable variant types in a single reference sample. With the input of recognized clinical genomics experts, the Seraseq ctDNA v4 Reference Materials focus on clinically relevant variants that provide actionable feedback about patient samples.

## FEATURES AND BENEFITS

- 93 unique multiplexed variants in 71 genes, covering 43 SNVs, 19 deletions, 5 insertions, 4 INDELS, 12 CNVs, and 10 translocations for the broadest coverage of clinically relevant variant types in a single reference material
- Four different variant allele frequencies (VAF) – WT (0%), 0.1%, 0.5%, and 5% - challenge your limit-of-detection or match typical AF of solid tumor NGS assays
- 83 variants are FDA drug targets
- Variants precisely quantitated with digital PCR and orthogonally validated by NGS
- Blended with well-characterized GM24385 human genomic DNA as background wild-type material
- Manufactured in cGMP-compliant, ISO 13485 certified facilities

## HIGHLIGHTS

Diverse set of 93 variants including CNVs and translocations

Improved manufacturing to reduce background noise

Consistently performing with batch-specific yield information provided

## ORDERING INFORMATION

Material #	Product	Concentration*	Fill Volume	Total Mass
0710-3097	Seraseq ctDNA Mutation Mix v4 AF0.1%	10 ng/μL	25 μL	250 ng
0710-3099	Seraseq ctDNA Mutation Mix v4 AF0.5%	10 ng/μL	25 μL	250 ng
0710-3100	Seraseq ctDNA Mutation Mix v4 AF5%	10 ng/μL	25 μL	250 ng
0710-3101	Seraseq ctDNA Mutation Mix v4 WT	10 ng/μL	25 μL	250 ng

\* Concentration targets are based on the Qubit dsDNA BR Assay.

## MUTATIONS PRESENT IN THE SERASEQ® ctDNA v4 PRODUCTS (0.1%, 0.5%, 5%)

Gene	Nucleotide Change	Transcript	Variant Type	COSMIC ID
AKT1	c.49G>A	NM_005163.2	SNV	COSM33765
AR	c.2623C>T	NM_000044.6	SNV	COSM238555
ATM	c.1058_1059del	NM_000051.4	Deletion	COSM21924
BRAF	c.1799T>A	NM_004333.6	SNV	COSM476
BRCA1	c.1961del	NM_007294.4	Deletion	COSM219054
BRCA2	c.7934del	NM_000059.4	Deletion	COSM1738241
CDKN2A	c.9_32dup	NM_000077.5	Insertion	COSM13442
CHEK1	c.676del	NM_001114122.3	Deletion	COSM1352376*
CHEK2	c.1116_1117delinsGT	NM_007194.4	INDEL	COSM384945
EGFR	c.2235_2249del	NM_005228.5	Deletion	COSM6223
EGFR	c.2303G>T	NM_005228.5	SNV	COSM6241
EGFR	c.2310_2311insGGT	NM_005228.5	Insertion	COSM12378
EGFR	c.2369C>T	NM_005228.5	SNV	COSM6240
EGFR	c.2389T>A	NM_005228.5	SNV	COSM6493937
EGFR	c.2573T>G	NM_005228.5	SNV	COSM6224
ERBB2	c.2313_2324dup	NM_004448.4	Insertion	COSM20959
ESR1	c.1613A>G	NM_000125.4	SNV	COSM94250
FGFR3	c.746C>G	NM_000142.5	SNV	COSM715
HRAS	c.182A>G	NM_005343.4	SNV	COSM499
HRAS	c.37G>C	NM_005343.4	SNV	COSM486
IDH1	c.394C>T	NM_005896.4	SNV	COSM28747
IDH2	c.419G>A	NM_002168.4	SNV	COSM41590
IDH2	c.515G>A	NM_002168.4	SNV	COSM33733
KIT	c.2361+67_2361+72delTTTTTT	NM_000222.3	Deletion	N/A
KIT	c.2447A>T	NM_000222.3	SNV	COSM1314
KRAS	c.183A>C	NM_004985.5	SNV	COSM554
KRAS	c.34G>T	NM_004985.5	SNV	COSM516
KRAS	c.35G>A	NM_004985.5	SNV	COSM521
MAP2K1	c.370C>T	NM_002755.4	SNV	COSM235614

Gene	Nucleotide Change	Transcript	Variant Type	COSMIC ID
MAP4K3	c.246-2475_246-2470delTTTTTT	NM_003618.4	NM_003618.4	N/A
MAP4K3	c.998-35_998-30delAAAAAA	NM_003618.4	Deletion	N/A
MET	c.3082+1del	NM_001127500.3	Deletion	COSM6947926
MLH1	c.232_243delinsATGTAAGG	NM_000249.4	INDEL	N/A
MSH2	c.1662-12_1677del	NM_000251.3	Deletion	N/A
MSH2	c.942+20_942+29delAAAAAAAAAA	NM_000251.3	Deletion	N/A
MSH6	c.2056_2060delinsCTTCTACCTCAAAA	NM_000179.3	INDEL	N/A
MTOR	c.6644C>A	NM_004958.4	SNV	COSM20417
NF1	c.3738_3747del	NM_001042492.3	Deletion	COSM510741
NRAS	c.182A>G	NM_002524.5	SNV	COSM584
NTRK1	c.1783G>A	NM_002529.4	SNV	COSM9113104
NTRK2	c.1915G>A	NM_006180.6	SNV	N/A
TERT	c.-124C>T	NM_198253.3	SNV	N/A
TERT	c.-146C>T	NM_198253.3	SNV	N/A
TP53	c.723del	NM_000546.6	Deletion	COSM6530
NTRK3	c.1867G>A	NM_001012338.3	SNV	COSM6951362
PALB2	c.839del	NM_024675.4	Deletion	COSM1376815
PDGFRA	c.2525A>T	NM_006206.6	SNV	COSM736
PIK3CA	c.1633G>A	NM_006218.4	SNV	COSM763
PIK3CA	c.3140A>G	NM_006218.4	SNV	COSM775
PIK3CA	c.3203dup	NM_006218.4	Insertion	COSM249879
PIK3R1	c.1727_1729del	NM_181523.3	Deletion	COSM35737
PMS2	c.861_864del	NM_000535.7	Deletion	COSM5547641
PTCH1	c.2307_2308delinsTT	NM_000264.5	INDEL	COSM17587
PTEN	c.800del	NM_000314.8	Deletion	COSM5809
PTEN	c.741dup	NM_000314.8	Insertion	COSM4986
RAD51C	c.242C>A	NM_058216.3	SNV	N/A
RAD51C	c.338dup	NM_058216.3	SNV	N/A
RAD51D	c.271A>T	NM_002878.4	SNV	N/A
RAD51D	c.392dup	NM_002878.4	SNV	N/A

Gene	Nucleotide Change	Transcript	Variant Type	COSMIC ID
RAF1	c.770C>T	NM_002880.4	SNV	COSM181063
RB1	c.751C>T	NM_000321.3	SNV	COSM878
RET	c.2753T>C	NM_020975.6	SNV	COSM965
SLC7A8	c.-231_-224delTTTTTTTT	NM_012244.4	Deletion	N/A
SMARCB1	c.118C>T	NM_003073.5	SNV	COSM1002
STK11	c.734+1G>T	NM_000455.5	SNV	COSM51523
RAD51D	c.392dup	NM_002878.4	SNV	N/A
RAF1	c.770C>T	NM_002880.4	SNV	COSM181063
RB1	c.751C>T	NM_000321.3	SNV	COSM878
RET	c.2753T>C	NM_020975.6	SNV	COSM965
SLC7A8	c.-231_-224delTTTTTTTT	NM_012244.4	Deletion	N/A
SMARCB1	c.118C>T	NM_003073.5	SNV	COSM1002
STK11	c.734+1G>T	NM_000455.5	SNV	COSM51523
TERT	c.-124C>T	NM_198253.3	SNV	N/A
TERT	c.-146C>T	NM_198253.3	SNV	N/A
TP53	c.723del	NM_000546.6	Deletion	COSM6530
TP53	c.743G>A	NM_000546.6	SNV	COSM10662
TP53	c.818G>A	NM_000546.6	SNV	COSM10660
TSC1	c.1263+1G>T	NM_000368.5	SNV	COSM1738312
TSC2	c.2640-1G>A	NM_000548.5	SNV	COSM3361675
VHL	c.481C>T	NM_000551.4	SNV	COSM17612
ZNF2	c.*1525_*1530delTTTTTT	NM_021088.4	Deletion	N/A
CD74::NRG1	Intron 6::Intron 5	NM_001025159.3::NM_013964.5	Translocation	N/A
CD74::ROS1	Intron 6::Intron 34	NM_001025159.3::NM_001378902.1	Translocation	N/A
COL1A1::PDGFB	Intron 25::Intron 1	NM_000088.3::NM_002608.3	Translocation	N/A
EML4::ALK	Intron 13::Intron 19	NM_019063.5::NM_004304.5	Translocation	N/A
ETV6::NTRK3	Intron 5::Intron 14	NM_001987.5::NM_002530.4	Translocation	N/A
FGFR2::BICC1	Intron 17::Intron 2	NM_000141.5::NM_001080512.3	Translocation	N/A
FGFR3::TACC3	Exon 18::Intron 7	NM_000142.5::NM_006342.3	Translocation	N/A
NCOA4::RET	Intron 7::Intron 11	NM_001145263.2::NM_020975.6	Translocation	N/A

Gene	Nucleotide Change	Transcript	Variant Type	COSMIC ID
PML::NTRK2	Intron 2::Intron 12	NM_002675.4::NM_006180.6	Translocation	N/A
TPM3::NTRK1	Intron 7::Intron 9	NM_153649.4::NM_002529.4	Translocation	N/A
AKT2	Amplification	19:40736224_40791252	CNV	N/A
CCND1	Amplification	NM_053056.3:11:69455924_69469242	CNV	N/A
CCNE1	Amplification	NM_001238.4:19:30302898_30315219	CNV	N/A
CDK4	Amplification	NM_000075.4:58141510_58146093	CNV	N/A
ERBB2	Amplification	NM_004448.4:17:37844347_37884911	CNV	N/A
FGF19	Amplification	NM_005117.3:11:69513006_69518790	CNV	N/A
FGF3	Amplification	NM_005247.4:11:69624736_69634184	CNV	N/A
FGF4	Amplification	NM_002007.4:69587797_69590109	CNV	N/A
FGFR1	Amplification	NM_023110.3:8:38268661_38326153	CNV	N/A
MET	Amplification	NM_001127500.3:7:116312250_116438431	CNV	N/A
MYC	Amplification	NM_002467.6:8:128747680_128755197	CNV	N/A
MYCN	Amplification	NM_005378.6:2:16080672_16087126	CNV	N/A

\* COSMIC uses transcript ENST00000427383.6

**NOTE:** Above list does not include variants present in the GM24385 background. Indels are defined as deletion/insertions less than 10 base pairs.



FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.  
 Seraseq® is a registered trademark of LGC Clinical Diagnostics, Inc. © 2024 LGC Clinical Diagnostics, Inc. All rights reserved.  
 MKT-01055-Rev. 01

seracare.com • 508.244.6400 • 800.676.1881 Toll Free • 508.634.3334 Fax