

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description:	SeraSeq® Myeloid Mutation DNA Mix		
Material Number:	0710-0408	Batch Number:	10724372
Material Description:	A multiplexed mixture of biosynthetic DNA targets with clinically relevant mutations precisely blended with a single, well-characterized genomic background at 15%, 10% or 5% allelic ratio		
Nominal Concentration:	15 ng/μL		
Avg Measured Concentration: (Qubit dsDNA BR kit)	20.2		
Fill Volume:	25 μL		
Date of Manufacture:	05-NOV-2024	Expiration Date:	05-NOV-2028
Test Methods:	Digital PCR (dPCR) using BioRad QX200™ Droplet Digital™ PCR System NGS using ArcherDx VariantPlex® Core Myeloid Kit, on ILMN MiSeq NGS System		

# Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Myeloid Mutation DNA Mix

Gene ID	HGVS	Protein	COSMIC ID	VAF Target (%)	VAF by dPCR (%)	VAF by NGS (%)
CALR	c.1099_1150del	p.L367fs*46	1738055	5	4.7	5.3
CSF3R	c.1853C>T	p.T618I	1737962	5	4.7	5.6
FLT3	c.1759_1800dup	p.N600_L601_insNEYVDFREYEYD	NA	5	4.9	2.8
IDH1	c.394C>T	p.R132C	28747	5	5.4	2.9
JAK2	c.1849G>T	p.V617F	12600	5	4.7	3.5
MPL	c.1544G>T	p.W515L	18918	5	4.9	4.6
NPM1	c.860_863dup	p.W288fs*12	17559	5	4.2	0.7
SF3B1	c.2098A>G	p.K700E	84677	5	4.1	4.6
SF3B1	c.1998G>T	p.K666N	131557	5	4.1	4.1
SRSF2	c.284_307del	p.P95_R102del	146289	5	4.7	8.3
ABL1	c.944C>T	p.T315I	12560	10	10.3	12.1
ASXL1	c.1900_1922del	p.E635fs*15	36165	10	8.2	8.1
ASXL1	c.1934dup	p.G646fs*12	34210	10	9.3	9.2
BRAF	c.1799T>A	p.V600E	476	10	9.7	8.2
CBL	c.1139T>C	p.L380P	34055	10	8.8	10.6
CBL	c.1259G>A	p.R420Q	34077	10	8.8	12.7
FLT3	c.1806_1807insGGGGCTTTCAGAGAATATG AATATGATCTCAAA)	K602_W603insGAFREYEYDLK	NA	10	8.6	0.7
FLT3	c.2503G>T	p.D835Y	783	10	10.3	8.7
JAK2	c.1624_1629del	p.N542_E543del	24440	10	13.1	8.8
MYD88	c.794T>C	p.L265P	85940	10	9.5	Not assayed
U2AF1	c.101C>T	p.S34F	166866	10	8.8	8.2
CEBPA	c.68dup	p.H24fs*84	18922	15	NA*	11.4
CEBPA	c.937_939dup	p.K313_V314insK	18099	15	NA*	10.7

\*Not assessed by dPCR, but confirmed by NGS (ArcherDx VariantPlex Core Myeloid)

Approval:

*Tur*

*11/26/2024*

Prepared By

Date