

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® CNV Breast Mix, +12 copies

Material Number: 0710-0413 Batch Number: 10713082

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Concentration (Qubit dsDNA BR Assay): Nominal value: 10 ng/μL
Average measured value: 15.5 ng/μL


Fill Volume: 20 μL

Date of Manufacture: 01 OCT 2024 Expiration Date: 01 OCT 2026

Test Method: BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System

Gene	Number of Additional Copies
ERBB2	12.7
FGFR3	11.5
MYC	12.5

Approval:



28 OCT 2024

Prepared By

Date