

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 18 Female - Matched Reference Material
Material Number: 0720-0172 Batch Number: 10724702
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
Concentration (Qubit dsDNA BR Assay): 33.5 ng/mL
Volume: 1.0 mL
Date of Manufacture: 09 DEC 2024 Expiration Date: 11 SEP 2027
Fetal Fraction by NGS (Chroma): 9.5%

Trisomy Call for Chromosome 18 by NGS (Chroma): Pass
Call for fetal gender by NGS (Chroma): XX
Storage: 2-8°C

Approval:

A handwritten signature in blue ink, appearing to be "MJA".

23-JAN-2025
Date